

ARG70451**Human MMP7 recombinant protein (Proenzyme) (His-tagged, C-ter)**

Package: 100 µg, 20 µg


Store at: -20°C

Summary

Product Description	E. coli expressed, His-tagged (C-ter) Human MMP7 recombinant protein
Tested Application	SDS-PAGE
Target Name	MMP7
Species	Human
A.A. Sequence	Met1 - Lys276
Expression System	E. coli
Alternate Names	MMP7; Matrix Metallopeptidase 7; Matrilysin; PUMP-1; MPSL1; Matrix Metalloproteinase 7 (Matrilysin, Uterine); Matrix Metalloproteinase-7; Uterine Metalloproteinase; Pump-1 Protease; Matrin; MMP-7

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images**ARG70451 Human MMP7 recombinant protein (Proenzyme) (His-tagged, C-ter) SDS-PAGE image**

SDS-PAGE analysis of ARG70451 Human MMP7 recombinant protein (Proenzyme) (His-tagged, C-ter)