

ARG70453 Human PGA5 recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Human PGA5 recombinant protein
Tested Application	SDS-PAGE
Target Name	PGA5
Species	Human
A.A. Sequence	Met1 - Ala388
Expression System	E. coli
Alternate Names	PGA5; Pepsinogen A5; Pepsinogen 5, Group I (Pepsinogen A); Pepsinogen-5; Pepsin A-5; Pepsin A; Pepsinogen-4

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 8.0)
Reconstitution	

Images



ARG70453 Human PGA5 recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70453 Human PGA5 recombinant protein (His-tagged, C-ter)