

ARG70454 Human OGG1 recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Human OGG1 recombinant protein
Tested Application	SDS-PAGE
Target Name	OGG1
Species	Human
A.A. Sequence	Met1 - Gly345
Expression System	E. coli
Alternate Names	OGG1; 8-Oxoguanine DNA Glycosylase; OGH1; MUTM; HOGG1; HMMH; 8-Hydroxyguanine DNA Glycosylase; N-Glycosylase/DNA Lyase; DNA-Apurinic Or Apyrimidinic Site Lyase

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70454 Human OGG1 recombinant protein (His-tagged, C-ter)
SDS-PAGE image

SDS-PAGE analysis of ARG70454 Human OGG1 recombinant protein (His-tagged, C-ter)