

ARG70456 Human Pleiotrophin recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Human PTN recombinant protein
Tested Application	SDS-PAGE
Target Name	PTN
Species	Human
A.A. Sequence	Gly33 - Asp168
Expression System	E. coli
Alternate Names	PTN; Pleiotrophin; HBNF; HBGF8; NEGF1; Heparin-Binding Neurite Outgrowth-Promoting Factor 1; Heparin-Binding Growth-Associated Molecule; Heparin-Binding Growth Factor 8; Heparin-Binding Brain Mitogen; Osteoblast-Specific Factor 1; HB-GAM; HBGF-8; HBNF-1; OSF-1; HBBM

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70456 Human Pleiotrophin recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70456 Human Pleiotrophin recombinant protein (His-tagged, C-ter)