

ARG70457 Mouse CDNF recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Mouse CDNF recombinant protein
Tested Application	SDS-PAGE
Target Name	CDNF
Species	Mouse
A.A. Sequence	Gln25 - Leu187
Expression System	E. coli
Alternate Names	CDNF; Cerebral Dopamine Neurotrophic Factor; Conserved Dopamine Neurotrophic Factor; ARMETL1; Arginine-Rich, Mutated In Early Stage Tumors-Like 1; ARMET-Like Protein 1

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70457 Mouse CDNF recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70457 Mouse CDNF recombinant protein (His-tagged, C-ter)