

ARG70459

Mouse Neurturin recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Mouse NRTN recombinant protein
Tested Application	SDS-PAGE
Target Name	NRTN
Species	Mouse
A.A. Sequence	Cys101 - Val195
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce proliferation in SH-SY5Y cells. The ED50 for this effect is < 50 ng/mL.
Alternate Names	NRTN; Neurturin; NTN; Prepro-Neurturin

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70459 Mouse Neurturin recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70459 Mouse Neurturin recombinant protein (Active) (His-tagged, N-ter)