

## ARG70460 Mouse Midkine recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (C-ter) Mouse Midkine recombinant protein
Tested Application	SDS-PAGE
Target Name	Midkine
Species	Mouse
A.A. Sequence	Gly31 - Val119
Expression System	E. coli
Alternate Names	MDK; Midkine; MK; NRG2; Neurite Outgrowth-Promoting Factor 2; Neurite Growth-Promoting Factor 2; Amphiregulin-Associated Protein; Midgestation And Kidney Protein

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70460 Mouse Midkine recombinant protein (His-tagged, C-ter)  
SDS-PAGE image

SDS-PAGE analysis of ARG70460 Mouse Midkine recombinant protein (His-tagged, C-ter)