

ARG70461

Human Galectin 2 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Human Galectin 2 recombinant protein
Tested Application	SDS-PAGE
Target Name	Galectin 2
Species	Human
A.A. Sequence	Thr2 - Glu132
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to agglutinate human red blood cells. The ED50 for this effect is < 20 µg/mL.
Alternate Names	LGALS2; Galectin 2; HL14; Lectin, Galactoside-Binding, Soluble, 2; Beta-Galactoside-Binding Lectin L-14-II; Lactose-Binding Lectin 2; S-Lac Lectin 2; Galectin-2; Gal-2

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70461 Human Galectin 2 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70461 Human Galectin 2 recombinant protein (Active) (His-tagged, N-ter)