

ARG70462

Human Galectin 8 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg

Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Human Galectin 8 recombinant protein
Tested Application	SDS-PAGE
Target Name	Galectin 8
Species	Human
A.A. Sequence	Met2 - Trp317
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to agglutinate human red blood cells. The ED50 for this effect is < 8 µg/mL.
Alternate Names	LGALS8; Galectin 8; PCTA-1; Lectin, Galactoside-Binding, Soluble, 8; Prostate Carcinoma Tumor Antigen 1 Po66 Carbohydrate-Binding Protein; Galectin-8; Po66-CBP; Gal-8

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70462 Human Galectin 8 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70462 Human Galectin 8 recombinant protein (Active) (His-tagged, N-ter)