

## ARG70463

## Human Galectin 9 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg

Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Human Galectin 9 recombinant protein
Tested Application	SDS-PAGE
Target Name	Galectin 9
Species	Human
A.A. Sequence	Ala2 - Thr355
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability of the immobilized protein to support the adhesion of Jurkat cells. The ED50 for this effect is < 3 µg/mL.
Alternate Names	LGALS9; Galectin 9; LGALS9A; Lectin, Galactoside-Binding, Soluble, 9; Tumor Antigen HOM-HD-21; Galectin-9; Ecalectin; Gal-9

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70463 Human Galectin 9 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70463 Human Galectin 9 recombinant protein (Active) (His-tagged, N-ter)