

## ARG70464 Human Galectin 10 recombinant protein (His-tagged, N-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

|                     |                                                                                                                                                                                                       |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Description | E. coli expressed, His-tagged (N-ter) Human Galectin 10 recombinant protein                                                                                                                           |
| Tested Application  | SDS-PAGE                                                                                                                                                                                              |
| Target Name         | Galectin 10                                                                                                                                                                                           |
| Species             | Human                                                                                                                                                                                                 |
| A.A. Sequence       | Ser2 - Arg142                                                                                                                                                                                         |
| Expression System   | E. coli                                                                                                                                                                                               |
| Alternate Names     | CLC; Charcot-Leyden Crystal Galectin; LGALS10 Gal-10; Charcot-Leyden Crystal Protein; Eosinophil Lysophospholipase; Lysolecithin Acylhydrolase; Lectin, Galactoside-Binding, Soluble, 10; Galectin-10 |

### Properties

|                   |                                                                               |
|-------------------|-------------------------------------------------------------------------------|
| Form              | Powder                                                                        |
| Purification Note | Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test. |
| Purity            | > 98% (by SDS-PAGE)                                                           |
| Buffer            | PBS (pH 7.4)                                                                  |
| Reconstitution    |                                                                               |

### Images



ARG70464 Human Galectin 10 recombinant protein (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70464 Human Galectin 10 recombinant protein (His-tagged, N-ter)