

ARG70465 Human Galectin 12 recombinant protein (His-tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Human Galectin 12 recombinant protein
Tested Application	SDS-PAGE
Target Name	Galectin 12
Species	Human
A.A. Sequence	Ser2 - Ser336
Expression System	E. coli
Alternate Names	LGALS12; Galectin 12; GRIP1; Galectin-Related Inhibitor Of Proliferation; Lectin, Galactoside-Binding, Soluble, 12; Galectin-12; Testicular Secretory Protein Li 26; Gal-12; GAL12

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70465 Human Galectin 12 recombinant protein (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70465 Human Galectin 12 recombinant protein (His-tagged, N-ter)