

ARG70466 Human Galectin 13 recombinant protein (His-tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Human Galectin 13 recombinant protein
Tested Application	SDS-PAGE
Target Name	LGALS13
Species	Human
A.A. Sequence	Ser2 - Asn139
Expression System	E. coli
Alternate Names	LGALS13; Galectin 13; PLAC8; PP13; Lectin, Galactoside-Binding, Soluble, 13; Galactoside-Binding Soluble Lectin 13; Placental Tissue Protein 13; Placental Protein 13; Gal-13; Beta-Galactoside-Binding Lectin; Galectin-13; GAL13

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70466 Human Galectin 13 recombinant protein (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70466 Human Galectin 13 recombinant protein (His-tagged, N-ter)