

ARG70467 Human Galectin 14 recombinant protein (His-tagged, N-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (N-ter) Human Galectin 14 recombinant protein
Tested Application	SDS-PAGE
Target Name	Galectin 14
Species	Human
A.A. Sequence	Ser2 - Asp139
Expression System	E. coli
Alternate Names	LGALS14; Galectin 14; PPL13; CLC2; Lectin, Galactoside-Binding, Soluble, 14; Charcot-Leyden Crystal Protein 2; Placental Protein 13-Like; Gal-14; Placental Protein 13-Like Protein; Galectin-14

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70467 Human Galectin 14 recombinant protein (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70467 Human Galectin 14 recombinant protein (His-tagged, N-ter)