

## ARG70471 Mouse CD70 recombinant protein (Active) (His-tagged, N-ter)

Package: 100 µg, 20 µg  
Store at: -20°C


### Summary

Product Description	E. coli expressed, His-tagged (N-ter) Active Mouse CD70 recombinant protein
Tested Application	SDS-PAGE
Target Name	CD70
Species	Mouse
A.A. Sequence	Gln47 - Pro195
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce proliferation in mouse T cells in the presence of the anti-CD3 antibody. The ED50 for this effect is < 4.5 µg/mL
Alternate Names	CD70; CD70 Molecule; CD27L; CD70 Antigen; CD27 Ligand; CD27-L; CD27LG; TNFSF7; Tumor Necrosis Factor Ligand Superfamily Member 7; Tumor Necrosis Factor (Ligand) Superfamily, Member 7

### Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

### Images



ARG70471 Mouse CD70 recombinant protein (Active) (His-tagged, N-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70471 Mouse CD70 recombinant protein (Active) (His-tagged, N-ter)