

ARG70472 Human Annexin A5 recombinant protein (His-tagged, C-ter)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, His-tagged (C-ter) Human Annexin A5 recombinant protein
Tested Application	SDS-PAGE
Target Name	Annexin A5
Species	Human
A.A. Sequence	Met1 - Asp320
Expression System	E. coli
Alternate Names	ANXA5; Annexin A5; CPB-I; PAP-I; Placental Anticoagulant Protein I; Vascular Anticoagulant-Alpha; Calphobindin I; Endonexin II; Annexin V; VAC-Alpha; RPRGL3; ANX5

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 98% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70472 Human Annexin A5 recombinant protein (His-tagged, C-ter) SDS-PAGE image

SDS-PAGE analysis of ARG70472 Human Annexin A5 recombinant protein (His-tagged, C-ter)