

ARG70473 Human IL18 recombinant protein (Active) (Tag free)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, Active Human IL18 recombinant protein
Tested Application	SDS-PAGE
Target Name	IL18
Species	Human
A.A. Sequence	Tyr37 - Asp193
Expression System	E. coli
Activity	Active
Activity Note	Determined by its ability to induce IFN gamma secretion in KG-1 cells. The ED50 for this effect is < 6 ng/mL.
Alternate Names	lbctadekin; Interferon gamma-inducing factor; Interleukin-18; IL-18; IFN-gamma-inducing factor; IL-1 gamma; IGIF; IL-1g; Interleukin-1 gamma; IL1F4

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70473 Human IL18 recombinant protein (Active) (Tag free) SDS-PAGE image

SDS-PAGE analysis of Human IL18 recombinant protein (Active) (Tag free)