

ARG70474 Human Activin A recombinant protein (Active) (Tag free)

Package: 100 µg, 20 µg
Store at: -20°C


Summary

Product Description	E. coli expressed, Active Human Activin A recombinant protein
Tested Application	SDS-PAGE
Target Name	Activin A
Species	Human
A.A. Sequence	Gly311-Ser426
Expression System	E. coli
Activity	Active
Alternate Names	Erythroid differentiation protein; Activin beta-A chain; FRP; EDF; Inhibin beta A chain

Properties

Form	Powder
Purification Note	Endotoxin level is < 0.1 EU/µg of the protein, as determined by the LAL test.
Purity	> 95% (by SDS-PAGE)
Buffer	PBS (pH 7.4)
Reconstitution	

Images



ARG70474 Human Activin A recombinant protein (Active) (Tag free)
SDS-PAGE image

SDS-PAGE analysis of ARG70474 Human Activin A recombinant protein (Active) (Tag free)