

ARG80208 Mouse Cxcl1 ELISA Kit

Package: 96 wells
Store at: 4°C

Component

Cat. No.	Component Name	Package	Temp
ARG80208-001	Antibody-coated microplate	8 X 12 strips	4°C. Unused strips should be sealed tightly in the air-tight pouch.
ARG80208-002	Standard (Lyophilized)	3 X 2 ng/vial	4°C
ARG80208-003	Standard diluent buffer	20 ml	4°C
ARG80208-004	Antibody conjugate concentrate	400 µl	4°C
ARG80208-005	Antibody diluent buffer	16 ml	4°C
ARG80208-006	HRP-Streptavidin concentrate	400 µl	4°C (Protect from light)
ARG80208-007	HRP-Streptavidin diluent buffer	16 ml	4°C
ARG80208-008	20X Wash buffer	50 ml	4°C
ARG80208-009	TMB substrate	12ml	4°C (Protect from light)
ARG80208-010	STOP solution	12ml	4°C
ARG80208-011	Plate sealer	4 strips	Room temperature

Summary

Product Description	ARG80208 Mouse Cxcl1 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Mouse Cxcl1 alpha in Serum, Plasma, Cell culture supernatants.
Tested Reactivity	Ms
Tested Application	ELISA
Specificity	No significant cross-reactivity or interference with Mouse MIP-1 alpha, MIP-1 beta, MIP-2, human GRO alpha, GRO beta, GRO gamma
Target Name	CXCL1
Species	Mouse
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm
Sensitivity	8 pg/ml
Sample Type	Serum, Plasma, Cell culture supernatants

Standard Range	15.6 - 1000 pg/ml
Sample Volume	100 µl
Precision	CV: less than 10%
Alternate Names	Growth-regulated alpha protein; SCYB1; Melanoma growth stimulatory activity; MGSA-a; GRO-alpha; GROa; NAP-3; FSP; C-X-C motif chemokine 1; Neutrophil-activating protein 3; GRO1; MGSA; IL-8 homologues (murine)

Application Instructions


Assay Time	4 hours
------------	---------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	Cxcl1
Gene Full Name	chemokine (C-X-C motif) ligand 1 (melanoma growth stimulating activity, alpha)
Background	This antimicrobial gene encodes a member of the CXC subfamily of chemokines. The encoded protein is a secreted growth factor that signals through the G-protein coupled receptor, CXC receptor 2. This protein plays a role in inflammation and as a chemoattractant for neutrophils. Aberrant expression of this protein is associated with the growth and progression of certain tumors. A naturally occurring processed form of this protein has increased chemotactic activity. Alternate splicing results in coding and non-coding variants of this gene. A pseudogene of this gene is found on chromosome 4. [provided by RefSeq, Sep 2014]
Function	Has chemotactic activity for neutrophils. May play a role in inflammation and exerts its effects on endothelial cells in an autocrine fashion. In vitro, the processed forms GRO-alpha(4-73), GRO-alpha(5-73) and GRO-alpha(6-73) show a 30-fold higher chemotactic activity. [UniProt]
Highlight	<p>Related products: Cxcl1 antibodies; Cxcl1 ELISA Kits;</p> <p>Related news: Immune signaling protein TLR4 has opposing roles in breast cancer development</p> <p>New ELISA data calculation tool: Simplify the ELISA analysis by GainData</p>
Research Area	Cancer kit; Cell Biology and Cellular Response kit; Immune System kit
PTM	N-terminal processed forms GRO-alpha(4-73), GRO-alpha(5-73) and GRO-alpha(6-73) are produced by proteolytic cleavage after secretion from peripheral blood monocytes.



ARG80208 Mouse Cxcl1 ELISA Kit standard curve image

ARG80208 Mouse Cxcl1 ELISA Kit results of a typical standard run with optical density reading at 450 nm.