

ARG80358 Human ANCA Screen ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG80358 Human ANCA Screen ELISA Kit is an ELISA test system for the determination of IgG class autoantibodies against PR3 and MPO in human serum and plasma (EDTA, heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	ANCA Screen
Conjugation	HRP
Sensitivity	diagn. 95.1 %
Sample Type	Serum and plasma (EDTA, heparin, citrate).
Standard Range	Cut - off
Sample Volume	10 µl

Application Instructions

Assay Time	30, 15, 15, 5 min (RT)
------------	------------------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: ANCA antibodies: ANCA ELISA Kits: New ELISA data calculation tool: Simplify the ELISA analysis by GainData
Research Area	Immune System kit