


**ARG80453**  
**Glutamate ELISA Kit**

Package: 96 wells  
Store at: 4°C

**Summary**

Product Description	ARG80453 Glutamate ELISA Kit is an Enzyme Immunoassay kit for the quantification of Glutamate in serum, plasma (EDTA) and urine.
Tested Reactivity	Hu, Ms, Rat
Tested Application	ELISA
Specificity	

**Images**



ARG80453 Human Glutamate ELISA Kit standard curve image

ARG80453 Human Glutamate ELISA Kit results of a typical standard.