

## ARG80455 Human Histamine ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG80455 Human Histamine ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human Histamine in plasma (EDTA, heparin) and urine.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Histamine
Conjugation	HRP
Sensitivity	0.12 ng/ml (plasma), 0.30 ng/ml (urine)
Sample Type	Plasma (EDTA, heparin) and urine.
Standard Range	0.5 - 50 ng/ml
Sample Volume	10 µl urine; 25 µl plasma (EDTA, heparin)

### Application Instructions


Assay Time	3 h, 30 min, 20 min (RT/shaker)
------------	---------------------------------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Highlight	Related products: <a href="#">Histamine antibodies: Histamine ELISA Kits:</a> New ELISA data calculation tool: <a href="#">Simplify the ELISA analysis by GainData</a>
Research Area	Immune System kit



ARG80455 Human Histamine ELISA Kit standard curve example image

The representative standard curve of ARG80455 Human Histamine ELISA Kit. The standard curve is for demonstration only and cannot be used in place of data generations at the time of assay. The standard curve should be generated each time the assay is performed.