


ARG80480 Human 5-HT / Serotonin ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG80480 Human 5-HT / Serotonin ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human 5-HT / Serotonin in serum, urine, platelets and platelet-rich plasma (PRP).
Tested Reactivity	Hu, Macaq, Mk
Tested Application	ELISA
Specificity	

Images



ARG80480 Human Serotonin ELISA Kit standard curve example image

The representative standard curve of ARG80480 Human Serotonin ELISA Kit. The standard curve is for demonstration only and cannot be used in place of data generations at the time of assay. The standard curve should be generated each time the assay is performed.