

ARG80538 Human Dengue Virus IgG antibody ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG80538 Human Dengue Virus IgG antibody ELISA Kit is an Enzyme Immunoassay kit for the determination of IgG class antibodies against Dengue virus in human serum and plasma (heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Specificity	The ELISA kit reacts to all serotypes of DENV-1, DENV-2, DENV-3, and DENV4
Target Name	Dengue Virus IgG antibody
Conjugation	HRP
Sensitivity	diagn. >90%
Sample Type	Serum and plasma (heparin, citrate).
Standard Range	Cut - off
Sample Volume	10 µl

Application Instructions

Assay Time	1 h (37°C), 30 min (RT), 15 min (RT/dark)
------------	---

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Dengue Virus antibodies ; Dengue Virus ELISA Kits ; Dengue Virus Duos / Panels ; Related news: Best NS1 antibodies in the market Fighting fire with fire: Genetically-engineered mosquitoes as an alternative to fight diseases Tools for studying Dengue Virus Exploring Antiviral Immune Response
Research Area	Immune System kit; Microbiology and Infectious Disease kit