

## ARG80545 Human EBV EBNA1 IgG antibody ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG80545 Human EBV EBNA1 IgG antibody ELISA Kit has been designed for the determination of specific IgG antibodies against EBV EBNA-1 in serum and plasma (EDTA, citrate and heparin).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	EBV EBNA1 IgG antibody
Conjugation	HRP
Sensitivity	0.62 U/ml
Sample Type	Serum and plasma (EDTA, citrate and heparin).
Standard Range	0.62 - 60 U/ml
Sample Volume	10 µl

### Application Instructions

Assay Time	1 h, 30, 15 min
------------	-----------------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Microbiology and Infectious Disease kit
---------------	---