

ARG80552 Human FSME (TBEV) IgG antibody ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG80552 Human FSME (TBEV) IgG antibody ELISA Kit is an Enzyme Immunoassay kit for the quantitative determination of FSME (TBEV) IgG antibody in human serum and plasma (citrate, heparin).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	FSME (TBEV) IgG antibody
Conjugation	HRP
Sensitivity	1.78 U/ml
Sample Type	Serum and plasma (EDTA, heparin).
Standard Range	50 - 300 U/ml
Sample Volume	10 µl

Application Instructions

Assay Time	1 h, 30, 20 min
------------	-----------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: FSME ELISA Kits: New ELISA data calculation tool: Simplify the ELISA analysis by GainData
Research Area	Immune System kit; Microbiology and Infectious Disease kit