

## ARG80579 Human Leishmania infantum IgG antibody ELISA Kit

Package: 96 wells

Store at: 4°C

### Summary

Product Description	ARG80579 Human Leishmania infantum IgG antibody ELISA Kit is an Enzyme Immunoassay kit for the determination of Leishmania infantum IgG antibody in serum and plasma (heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Leishmania infantum IgG antibody
Conjugation	HRP
Sensitivity	diagn. 91 %.
Sample Type	Serum and plasma (heparin, citrate).
Standard Range	Cut - off
Sample Volume	10 µl serum or plasma

### Application Instructions

Assay Time	1 h (37°C),30, 15 min (RT)
------------	----------------------------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Immune System kit
---------------	-------------------