

## ARG80597 Human Parvovirus B19 IgM antibody ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG80597 Human Parvovirus B19 IgM antibody ELISA Kit provides materials for the determination of IgM-class antibodies to Parvovirus B19 in serum and plasma (EDTA, heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Parvovirus B19 IgM antibody
Conjugation	HRP
Sensitivity	diagn. 100%
Sample Type	Serum and plasma (EDTA, heparin, citrate).
Standard Range	cut - off standard
Sample Volume	100 µl

### Application Instructions

Assay Time	1 h, 30, 15 min
------------	-----------------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Microbiology and Infectious Disease kit
---------------	-----------------------------------------