

ARG80611 Human Toxocara canis IgG antibody ELISA Kit

Package: 96 wells

Store at: 4°C

Summary

Product Description	ARG80611 Human Toxocara canis IgG antibody ELISA Kit is intended for the determination of antibodies against Toxocara canis circulating in human serum and plasma (heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Toxocara canis IgG antibody
Conjugation	HRP
Sensitivity	> 98%
Sample Type	Serum and plasma (heparin, citrate).
Standard Range	cut - off index
Sample Volume	10 µl

Application Instructions

Assay Time	1 h, 30, 15 min
------------	-----------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Immune System kit; Microbiology and Infectious Disease kit
---------------	--