

ARG80652 Mouse/Rat Corticosterone ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG80652 Mouse/Rat Corticosterone ELISA Kit is a competitive immunoassay for the quantification of corticosterone in rat and mouse serum, plasma and cell culture medium.
Tested Reactivity	Ms, Rat
Tested Application	ELISA
Target Name	Corticosterone
Conjugation	HRP
Sensitivity	6.1 ng/ml
Sample Type	Serum, plasma and cell culture medium.
Standard Range	15 - 2250 ng/ml
Sample Volume	10 µl
Precision	Intra-Assay CV: 7.4% Inter-Assay CV: 7.6%

Application Instructions


Assay Time	2 h (shaking), 30 min
------------	-----------------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Corticosterone ELISA Kits; New ELISA data calculation tool: Simplify the ELISA analysis by GainData
Research Area	Signaling Transduction kit



ARG80652 Mouse/Rat Corticosterone ELISA Kit standard curve image

ARG80652 Mouse/Rat Corticosterone ELISA Kit results of a typical standard.