

ARG80756 Human NSE ELISA Kit

Package: 96 wells
Store at: 4°C

Summary


Product Description	ARG80756 Human NSE ELISA Kit is intended for the quantification of NSE in human serum.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	NSE
Conjugation	HRP
Sensitivity	0.19 ng/ml
Sample Type	Serum.
Standard Range	7 - 135 ng/ml
Sample Volume	25 µl
Alternate Names	Neural enolase; NSE; Enolase 2; Gamma-enolase; 2-phospho-D-glycerate hydro-lyase; HEL-S-279; Neuron-specific enolase; EC 4.2.1.11

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	ENO2
Gene Full Name	enolase 2 (gamma, neuronal)
Background	



ARG80756 Human NSE ELISA Kit standard curve image

ARG80756 Human NSE ELISA Kit results of a typical standard run with optical density reading at 450 nm.