


ARG81240 proBDNF ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG81240 proBDNF ELISA Kit is an Enzyme Immunoassay kit for the quantification of proBDNF in Human and rat cell culture supernatants and cell lysates; human serum and plasma (EDTA, citrate) and Rat tissue extracts.
Tested Reactivity	Hu, Rat
Predict Reactivity	Ms
Tested Application	ELISA
Specificity	

Images



ARG81240 proBDNF ELISA Kit standard curve image

ARG81240 proBDNF ELISA Kit results of a typical standard run with optical density reading at 450 nm.