

ARG81360 Human CD279 / PD-1 ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG81360 Human CD279 / PD-1 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human CD279 / PD-1 in serum, plasma and cell culture supernatant.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	CD279 / PD-1
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	0.44 pg/ml
Sample Type	Serum, plasma and cell culture supernatant.
Standard Range	15.6 - 1000 pg/ml
Sample Volume	100 µl
Alternate Names	hPD-I; CD279; PD-1; Protein PD-1; CD antigen CD279; PD1; hSLE1; SLEB2; Programmed cell death protein 1; hPD-1

Application Instructions

Assay Time	~ 3.5 hours
------------	-------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	PDCD1
Gene Full Name	programmed cell death 1
Background	CD279 / PD-1 is a cell surface membrane protein of the immunoglobulin superfamily. This protein is expressed in pro-B-cells and is thought to play a role in their differentiation. In mice, expression of this gene is induced in the thymus when anti-CD3 antibodies are injected and large numbers of thymocytes undergo apoptosis. Mice deficient for this gene bred on a BALB/c background developed dilated cardiomyopathy and died from congestive heart failure. These studies suggest that this gene product may also be important in T cell function and contribute to the prevention of autoimmune diseases. [provided by RefSeq, Jul 2008]

Function

CD279 / PD-1 is an inhibitory receptor on antigen activated T-cells. It plays a critical role in induction and maintenance of immune tolerance to self (PubMed:21276005). Delivers inhibitory signals upon binding to ligands CD274/PDCD1L1 and CD273/PDCD1LG2 (PubMed:21276005). Following T-cell receptor (TCR) engagement, PDCD1 associates with CD3-TCR in the immunological synapse and directly inhibits T-cell activation. Suppresses T-cell activation through the recruitment of PTPN11/SHP-2: following ligand-binding, PDCD1 is phosphorylated within the ITSM motif, leading to the recruitment of the protein tyrosine phosphatase PTPN11/SHP-2 that mediates dephosphorylation of key TCR proximal signaling molecules, such as ZAP70, PRKQC/PKTheta and CD247/CD3zeta.

The PDCD1-mediated inhibitory pathway is exploited by tumors to attenuate anti-tumor immunity and escape destruction by the immune system, thereby facilitating tumor survival (PubMed:28951311). The interaction with CD274/PDCD1L1 inhibits cytotoxic T lymphocytes (CTLs) effector function (PubMed:28951311). The blockage of the PDCD1-mediated pathway results in the reversal of the exhausted T-cell phenotype and the normalization of the anti-tumor response, providing a rationale for cancer immunotherapy (PubMed:22658127, PubMed:25034862, PubMed:25399552). [UniProt]

Highlight

Related products:

[PD-1 antibodies](#); [PD-1 ELISA Kits](#); [PD-1 Duos / Panels](#);

Related news:

[Detecting exosomal PD-L1 secreted by cancer cells](#)


[The best solution for PD-1/PD-L1 research](#)

[Examining CTL/NK-mediated cytotoxicity by ELISA](#)

New ELISA data calculation tool:


[Simplify the ELISA analysis by GainData](#)

Images



ARG81360 Human CD279 / PD-1 ELISA Kit standard curve image

ARG81360 Human CD279 / PD-1 ELISA Kit results of a typical standard run with optical density reading at 450 nm. (using CALIBRATOR DILUENT I)



ARG81360 Human CD279 / PD-1 ELISA Kit standard curve image

ARG81360 Human CD279 / PD-1 ELISA Kit results of a typical standard run with optical density reading at 450 nm. (using CALIBRATOR DILUENT II)