

ARG81664 Mouse Neurotrophin 3 ELISA Kit


Package: 96 wells
Store at: 4°C

Component

Cat. No.	Component Name	Package	Temp
ARG81664-001	Antibody-coated microplate	8 X 12 strips	4°C. Unused strips should be sealed tightly in the air-tight pouch.
ARG81664-002	Standard	2 X 10 ng/vial	4°C
ARG81664-003	Standard/Sample diluent	30 ml (Ready to use)	4°C
ARG81664-004	Antibody conjugate concentrate (100X)	1 vial (100 µl)	4°C
ARG81664-005	Antibody diluent buffer	12 ml (Ready to use)	4°C
ARG81664-006	HRP-Streptavidin concentrate (100X)	1 vial (100 µl)	4°C
ARG81664-007	HRP-Streptavidin diluent buffer	12 ml (Ready to use)	4°C
ARG81664-008	25X Wash buffer	20 ml	4°C
ARG81664-009	TMB substrate	10 ml (Ready to use)	4°C (Protect from light)
ARG81664-010	STOP solution	10 ml (Ready to use)	4°C
ARG81664-011	Plate sealer	4 strips	Room temperature

Summary

Product Description	ARG81664 Mouse Neurotrophin 3 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Mouse Neurotrophin 3 in serum and cell culture supernatants.
Tested Reactivity	Ms
Predict Reactivity	Bov, Chk, Dog, Hrs, Mk, Rb
Tested Application	ELISA
Specificity	



ARG81664 Mouse Neurotrophin 3 ELISA Kit standard curve image

ARG81664 Mouse Neurotrophin 3 ELISA Kit results of a typical standard run with optical density reading at 450 nm.