

ARG82015 Polyphenol Oxidase Assay Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG82015 Polyphenol Oxidase Assay Kit is a detection kit for the quantification of Polyphenol Oxidase Activity in urine, serum, plasma, tissue extracts, cell lysate and cell culture supernatants.
Tested Reactivity	All
Tested Application	FuncSt
Target Name	Polyphenol Oxidase
Conjugation Note	Read at 410 nm.
Sample Type	Urine, serum, plasma, tissue extracts, cell lysate and cell culture supernatants.
Sample Volume	50 µl

Application Instructions

Assay Time	~ 10 min
------------	----------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.