

ARG82137
Antioxidant Assay Kit (Colorimetric)Package: 100 tests
Store at: -20°C

Summary


Product Description	ARG82137 Antioxidant Assay Kit (Colorimetric) is a detection kit for the quantification of Antioxidant in serum, plasma, urine, saliva, food and beverage.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Antioxidant
Conjugation Note	Read at 570 nm
Sensitivity	1.5 µM
Detection Range	1.5 - 1000 µM Trolox equivalents
Sample Type	serum, plasma, urine, saliva, food and beverage
Sample Volume	20 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	20 min

Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82137 Antioxidant Assay Kit (Colorimetric) typical data demonstration image

ARG82137 Antioxidant Assay Kit (Colorimetric) results of a typical data with optical density reading at 570 nm.