

## ARG82177 Magnesium Assay Kit (Colorimetric)

Package: 250 tests  
Store at: 4°C

### Summary


Product Description	ARG82177 Magnesium Assay Kit (Colorimetric) is a detection kit for the quantification of Magnesium in Biological, food and beverage.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Magnesium
Conjugation Note	Read at 500 nm
Sensitivity	0.1 mg/dL (41 µM)
Detection Range	0.1 mg/dL (41 µM) - 3 mg/dL (1.2 mM)
Sample Type	Biological, food and beverage
Sample Volume	5 µl

### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	10 min

### Properties

Form	Liquid
Storage instruction	Store the kit at 2-8°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82177 Magnesium Assay Kit (Colorimetric) typical data demonstration image

ARG82177 Magnesium Assay Kit (Colorimetric) results of a typical data with optical density reading at 500 nm.