

**ARG82206**  
**GFP ELISA Kit**Package: 96 wells  
Store at: 4°C, -80°C

### Summary

Product Description	ARG82206 GFP ELISA Kit is an Enzyme Immunoassay kit for the quantification of Green fluorescent protein (GFP) or GFP fusion protein in cell culture supernatants, cell lysate and tissue lysate.
Tested Reactivity	Other
Tested Application	ELISA
Specificity	Detect GFP, BFP, CFP, and YFP from <i>Aequorea victoria</i> .
Target Name	GFP
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	30 pg/ml
Sample Type	Cell culture supernatants, cell lysate and tissue lysate.
Standard Range	31.2 - 2000 pg/ml
Sample Volume	100 µl

### Application Instructions


Assay Time	5.5 hours
------------	-----------

### Properties

Form	96 well
Storage instruction	Store components at 4°C or -80°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Highlight	Related products: <a href="#">GFP antibodies</a> ; <a href="#">GFP ELISA Kits</a> ; <a href="#">GFP Duos / Panels</a> ; New ELISA data calculation tool: <a href="#">Simplify the ELISA analysis by GainData</a>
-----------	---



ARG82206 GFP ELISA Kit standard curve image

ARG82206 GFP ELISA Kit results of a typical standard run with optical density reading at 450 nm.

---