

**ARG82226**  
**Bile Acid Assay Kit (Colorimetric)**Package: 100 assay  
Store at: -20°C, -80°C

### Summary


Product Description	ARG82226 Bile Acid Assay Kit (Colorimetric) is a detection kit for the quantification of Bile Acid in serum, plasma, cell lysate and tissue lysate.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Bile Acid
Conjugation Note	Read at 405 nm.
Sensitivity	1 µM
Sample Type	Serum, plasma, cell lysate and tissue lysate.
Standard Range	0.78 - 25 µM
Sample Volume	20 µl

### Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	~ 1 hour

### Properties

Form	Liquid
Storage instruction	Store components at -20°C or -80°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82226 Bile Acid Assay Kit (Colorimetric) standard curve image

ARG82226 Bile Acid Assay Kit (Colorimetric) results of a typical standard run with optical density reading at 405 - 630 nm.