

ARG82245 Human BNP ELISA Kit


Package: 96 wells

Store at: 4°C

Summary

Product Description	ARG82245 Human BNP ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human BNP in serum and plasma. (Peptide extraction is recommended)
Tested Reactivity	Hu
Tested Application	ELISA
Specificity	

Images



ARG82245 Human BNP ELISA Kit standard curve image

ARG82245 Human BNP ELISA Kit results of a typical standard run with optical density reading at 450 nm.