

ARG82271 Rat IgG2b ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG82271 Rat IgG2b ELISA Kit is an Enzyme Immunoassay kit for the quantification of Rat IgG2b in serum, plasma and cell culture supernatants.
Tested Reactivity	Rat
Tested Application	ELISA
Target Name	IgG2b
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	0.4 ng/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	0.78 - 50 ng/ml
Sample Volume	100 µl
Precision	Intra-Assay CV: less than 10% Inter-Assay CV: less than 10%

Application Instructions


Assay Time	~ 3.5 hours
------------	-------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: IgG2 ELISA Kits; New ELISA data calculation tool: Simplify the ELISA analysis by GainData
-----------	--



ARG82271 Rat IgG2b ELISA Kit standard curve image

ARG82271 Rat IgG2b ELISA Kit results of a typical standard run with optical density reading at 450 nm.