

ARG82284 Mouse MMP13 ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG82284 Mouse MMP13 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Mouse MMP13 in serum, plasma and cell culture supernatant.
Tested Reactivity	Ms
Tested Application	ELISA
Specificity	Not react with following recombinant mouse proteins: Adiponectin, ApoAI, BMP1, BMP2, BMP3, BMP4, BMP5, BMP6, BMP7, CCL2, CCL4, CCL5, CRP, FGF acidic, HGF, HSP27, IGF1, IL1 alpha, IL1 beta, IL1RA, IL4, IL5, IL6, IL8, IL10, IL12, IL13, IL15, IL17C, MMP131, MMP133, IFN alpha, IFN beta, IFN gamma, MMP2, MMP3, MMP9, PDGF, PLA2G7, prolactin, TLR1, TLR2, TLR3, TLR4, TLR9, TGF beta 1, TGF beta 2, TGF beta 3, TNF alpha, TNF RI, TNF RII, VEGF, VEGFR1.
Target Name	MMP13
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	31 pg/ml
Sample Type	Serum, plasma and cell culture supernatant.
Standard Range	62 - 4000 pg/ml
Sample Volume	100 µl
Precision	Intra-Assay CV: 5% Inter-Assay CV: 10%
Alternate Names	CLG3; EC 3.4.24.-; MANDP1; MMP-13; Collagenase 3; Matrix metalloproteinase-13

Application Instructions

Assay Time	~ 3 hours
------------	-----------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	MMP13
Gene Full Name	matrix metalloproteinase 13
Background	Proteins of the matrix metalloproteinase (MMP) family are involved in the breakdown of extracellular

matrix in normal physiological processes, such as embryonic development, reproduction, and tissue remodeling, as well as in disease processes, such as arthritis and metastasis. Most MMP's are secreted as inactive propeptides which are activated when cleaved by extracellular proteinases. The protein encoded by this gene cleaves type II collagen more efficiently than types I and III. It may be involved in articular cartilage turnover and cartilage pathophysiology associated with osteoarthritis. The gene is part of a cluster of MMP genes which localize to chromosome 11q22.3. [provided by RefSeq, Jul 2008]

Function

Plays a role in the degradation of extracellular matrix proteins including fibrillar collagen, fibronectin, TNC and ACAN. Cleaves triple helical collagens, including type I, type II and type III collagen, but has the highest activity with soluble type II collagen. Can also degrade collagen type IV, type XIV and type X. May also function by activating or degrading key regulatory proteins, such as TGFB1 and CTGF. Plays a role in wound healing, tissue remodeling, cartilage degradation, bone development, bone mineralization and ossification. Required for normal embryonic bone development and ossification. Plays a role in the healing of bone fractures via endochondral ossification. Plays a role in wound healing, probably by a mechanism that involves proteolytic activation of TGFB1 and degradation of CTGF. Plays a role in keratinocyte migration during wound healing. May play a role in cell migration and in tumor cell invasion. [UniProt]

Highlight

Related products:
[MMP13 antibodies](#); [MMP13 ELISA Kits](#);
Related news:
[Detecting MMPs and their non-ECM substrates](#)
New ELISA data calculation tool:
[Simplify the ELISA analysis by GainData](#)

PTM

The proenzyme is activated by removal of the propeptide; this cleavage can be effected by other matrix metalloproteinases, such as MMP2, MMP3 and MMP14 and may involve several cleavage steps. Cleavage can also be autocatalytic, after partial maturation by another protease or after treatment with 4-aminophenylmercuric acetate (APMA) (in vitro).


N-glycosylated.

Tyrosine phosphorylated by PKDCC/VLK. [UniProt]

Cellular Localization

Secreted, extracellular space, extracellular matrix. Secreted. [UniProt]

Images



ARG82284 Mouse MMP13 ELISA Kit standard curve image

ARG82284 Mouse MMP13 ELISA Kit results of a typical standard run with optical density reading at 450 nm.