

ARG82298 Rat IgM ELISA Kit

Package: 96 wells
Store at: 4°C

Summary


Product Description	ARG82298 Rat IgM ELISA Kit is an Enzyme Immunoassay kit for the quantification of Rat IgM in serum, plasma and cell culture supernatants.
Tested Reactivity	Rat
Tested Application	ELISA
Target Name	IgM
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	0.16 ng/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	0.313 - 20 ng/ml
Sample Volume	100 µl
Precision	CV: less than 8.5%

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: IgM antibodies: IgM ELISA Kits: IgM Duos / Panels: New ELISA data calculation tool: Simplify the ELISA analysis by GainData
-----------	--



ARG82298 Rat IgM ELISA Kit standard curve image

ARG82298 Rat IgM ELISA Kit results of a typical standard run with optical density reading at 450 nm.
