

ARG82508 Mouse CNTN6 / Contactin 6 ELISA Kit

Package: 96 wells
Store at: 4°C

Component

Cat. No.	Component Name	Package	Temp
ARG82508-001	Antibody-coated microplate	8 X 12 strips	4°C. Unused strips should be sealed tightly in the air-tight pouch.
ARG82508-002	Standard	2 X 10 ng/vial	4°C
ARG82508-003	Standard/Sample diluent	30 ml (Ready to use)	4°C
ARG82508-004	Antibody conjugate concentrate (100X)	1 vial (100 µl)	4°C
ARG82508-005	Antibody diluent buffer	12 ml (Ready to use)	4°C
ARG82508-006	HRP-Streptavidin concentrate (100X)	1 vial (100 µl)	4°C
ARG82508-007	HRP-Streptavidin diluent buffer	12 ml (Ready to use)	4°C
ARG82508-008	25X Wash buffer	20 ml	4°C
ARG82508-009	TMB substrate	10 ml (Ready to use)	4°C (Protect from light)
ARG82508-010	STOP solution	10 ml (Ready to use)	4°C
ARG82508-011	Plate sealer	4 strips	Room temperature

Summary

Product Description	ARG82508 Mouse CNTN6 / Contactin 6 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Mouse CNTN6 / Contactin 6 in serum, plasma (EDTA, heparin, citrate) and cell culture supernatants.
Tested Reactivity	Ms
Tested Application	ELISA
Target Name	CNTN6 / Contactin 6
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	78 pg/ml
Sample Type	Serum, plasma (EDTA, heparin, citrate) and cell culture supernatants.
Standard Range	156 - 10000 pg/ml
Sample Volume	100 µl
Precision	Intra-Assay CV: 4.4% Inter-Assay CV: 6.7%

Alternate Names NB3; hNB-3; Neural recognition molecule NB-3; Contactin-6

Application Instructions

Assay Time ~ 5 hours


Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	CNTN6
Gene Full Name	contactin 6
Background	The protein encoded by this gene is a member of the immunoglobulin superfamily. It is a glycosylphosphatidylinositol (GPI)-anchored neuronal membrane protein that functions as a cell adhesion molecule. It may play a role in the formation of axon connections in the developing nervous system. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]
Function	Contactins mediate cell surface interactions during nervous system development. Participates in oligodendrocytes generation by acting as a ligand of NOTCH1. Its association with NOTCH1 promotes NOTCH1 activation through the released notch intracellular domain (NICD) and subsequent translocation to the nucleus. Involved in motor coordination (By similarity). [UniProt]
Highlight	Related products: Contactin antibodies ; Contactin ELISA Kits ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData
Cellular Localization	Cell membrane; Lipid-anchor, GPI-anchor. [UniProt]

Images



ARG82508 Mouse CNTN6 / Contactin 6 ELISA Kit standard curve image

ARG82508 Mouse CNTN6 / Contactin 6 ELISA Kit results of a typical standard run with optical density reading at 450 nm.