

ARG82670 Monkey IFN alpha (pan) ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG82670 Monkey IFN alpha (pan) ELISA Kit is an Enzyme Immunoassay kit for the quantification of Monkey IFN alpha (pan) in serum, plasma and cell culture supernatants.
Tested Reactivity	Mk
Tested Application	ELISA
Target Name	IFN alpha (pan)
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	5.5 pg/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	10.9 - 700 pg/ml
Sample Volume	100 µl
Precision	CV: less than 10%
Alternate Names	IFNA1; Interferon Alpha 1; IFN-AlphaD; IFN-ALPHA IFNA13; IFNA@; IFL; IFN; Interferon Alpha-1/13; Interferon Alpha 1b; Interferon Alpha-D; IFN-Alpha-1/13; LeIF D; Interferon-Alpha1; IFN-Alpha 1b

Application Instructions

Assay Time	~ 3.5 hours
------------	-------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	IFNA1
Gene Full Name	interferon, alpha 1
Background	The protein encoded by this gene is produced by macrophages and has antiviral activity. This gene is intronless and the encoded protein is secreted. [provided by RefSeq, Sep 2011]
Function	Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase. [UniProt]
Highlight	Related products:

[IFN alpha antibodies: IFN alpha ELISA Kits:](#)

Related news:

[Exploring Antiviral Immune Response](#)

[circNDUFB2, a circular RNA \(circRNA\), activates anti-tumor immunity](#)


New ELISA data calculation tool:

[Simplify the ELISA analysis by GainData](#)

Cellular Localization

Secreted. [UniProt]

Images



ARG82670 Monkey IFN alpha (pan) ELISA Kit standard curve image

ARG82670 Monkey IFN alpha (pan) ELISA Kit results of a typical standard run with optical density reading at 450 nm.