


ARG82789 Human Metanephrine ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG82789 Human Metanephrine ELISA Kit is an Enzyme Immunoassay kit for the quantification of Metanephrine in Human plasma (EDTA, heparin).
Tested Reactivity	Hu
Tested Application	ELISA
Specificity	

Images



ARG82789 Human Metanephrine ELISA Kit standard curve image

ARG82789 Human Metanephrine ELISA Kit results of a typical standard run with optical density reading at 450 nm.