

ARG82794 Human 2-MET (Metanephrine + Normetanephrine) ELISA Kit


Package: 2x96 wells

Store at: 4°C

Summary


Product Description	ARG82794 Human 2-MET (Metanephrine + Normetanephrine) ELISA Kit is an Enzyme Immunoassay kit for the quantification of Metanephrine and Normetanephrine in Human urine.
Tested Reactivity	Hu
Tested Application	ELISA
Specificity	

Images



ARG82794 Human 2-MET (Metanephrine + Normetanephrine) ELISA Kit standard curve image

ARG82794 Human 2-MET (Metanephrine + Normetanephrine) ELISA Kit results of a typical standard run with optical density reading at 450 nm.



ARG82794 Human 2-MET (Metanephrine + Normetanephrine) ELISA Kit standard curve image

ARG82794 Human 2-MET (Metanephrine + Normetanephrine) ELISA Kit results of a typical standard run with optical density reading at 450 nm.