

**ARG82822**  
**HIV1 p24 ELISA Kit**Package: 96 wells  
Store at: 4°C, -20°C

### Summary


Product Description	ARG82822 HIV1 p24 ELISA Kit is an Enzyme Immunoassay kit for the quantification of HIV1 p24 in lentiviral cell culture supernatant.
Tested Reactivity	HIV
Tested Application	ELISA
Target Name	HIV1 p24
Conjugation	HRP
Conjugation Note	Read at 450 nm.
Sensitivity	1 ng/ml
Sample Type	Lentiviral cell culture supernatant.
Standard Range	1.56 - 100 ng/ml
Sample Volume	110 µl

### Application Instructions

Assay Time	Overnight, and 2.5 hours.
------------	---------------------------

### Properties

Form	96 well
Storage instruction	Store components at 4°C or -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82822 HIV1 p24 ELISA Kit standard curve image

ARG82822 HIV1 p24 ELISA Kit results of a typical standard run with optical density reading at 450 nm.