

ARG82844 Mouse IgG3 ELISA kit

Package: 96 wells
Store at: 4°C


Summary

Product Description	ARG82844 Mouse IgG3 ELISA kit is an Enzyme Immunoassay kit for the quantification of Mouse IgG3 in serum, plasma and cell culture supernatants.
Tested Reactivity	Ms
Tested Application	ELISA
Target Name	IgG3
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	0.78 ng/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	1.56 - 100 ng/ml
Sample Volume	100 µl

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Images



ARG82844 Mouse IgG3 ELISA kit standard curve image

ARG82844 Mouse IgG3 ELISA kit results of a typical standard run with optical density reading at 450 nm