

## ARG82908 Human Rubella IgG antibody ELISA Kit

Package: 96 wells

Store at: 4°C

### Summary

Product Description	ARG82908 Human Rubella IgG antibody ELISA Kit has been designed for the quantification of specific IgG antibodies against Rubella in serum and plasma.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Rubella IgG antibody
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 451 nm
Sensitivity	0.45 IU/mL
Sample Type	Serum and plasma (Citrate, heparin).
Standard Range	10 - 100 IU/mL
Sample Volume	5-10 µl

### Application Instructions

Assay Time	1 h, 30, 15 min
------------	-----------------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	for laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Highlight	Related products: <a href="#">Rubella antibodies</a> ; <a href="#">Rubella ELISA Kits</a> ; New ELISA data calculation tool: <a href="#">Simplify the ELISA analysis by GainData</a>
Research Area	Microbiology and Infectious Disease kit